2021 2nd

OSSTEM IMPLANT

**CONSENSUS  
REPORT**

## Contents

### Preface

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
Reference

Presented by Dr. 김기성

• KJ Lee, et. al. Influence of crown-to-implant ratio on peri-implant marginal bone loss in the posterior region: a five-year retrospective study. J Periodontal Implant Sci. 2012;42:231-236

• Yiman Tang, et. al. Influence of crown-to-implant ratio and different prosthetic designs on the clinical conditions of short implants in posterior regions: A 4-year retrospective clinical and radiographic study. Clin Implant Dent Relat Res. 2020;22:119-127

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## Reference

- 보건복지부, et. al. 치과감염관리 표준정책 매뉴얼. 2020
- Formlabs, Instructions for Use SURGICAL GUIDE RESIN, 11/12/2019 PRNT-0001, Rev: 01
- Gréta Török, et. al. Effects of disinfection and sterilization on the dimensional changes and mechanical properties of 3D printed surgical guides for implant therapy – pilot study. BMC Oral Health. (2020) 20:19
- Peter N Smith, et. al. Microbial contamination and the sterilization/disinfection of surgical guides used in the placement of endosteal implants. Int J Oral Maxillofac Implants. Mar-Apr 2011;26(2):274-81

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천공 (Perforation)

## Reference

- Liat Chaushu, et. al. Sinus Augmentation with Simultaneous, Non-Submerged, Implant Placement Using a Minimally Invasive Hydraulic Technique. Medicina. 2020, 56, 75

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
#### Reference

- Tiziano Testori, et. al. Perforation Risk Assessment in Maxillary Sinus Augmentation with Lateral Wall Technique. The International Journal of Periodontics & Restorative Dentistry. Volume 40, Number 3, 2020
- SC Wen, et. al. The influence of sinus membrane thickness upon membrane perforation during transcrestal sinus lift procedure. Clin Oral Implants Res. 2015 : 26(10) : 1158-1164
- YH Lin, et. al. The influence of sinus membrane thickness upon membrane perforation during transcrestal sinus lift procedure. Clin Oral Implants Res. 2016 : 27(5) : 612-617
- Alberto Monje, et. al. Schneiderian Membrane Thickness and Clinical Implications for Sinus Augmentation: A Systematic Review and Meta-Regression Analyses. J Periodontol. 2016 : 87(8) : 888-899
- D Uçar Boyacıgil, et. al. The effect of residual bone height and membrane thickness on sinus membrane perforation in crestal sinus grafting: A prospective clinical study. International Journal of Oral and Maxillofacial Surgery. 2021 : 50(2) : 251-257

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
### Reference

- Fouad Khoury, et. al. Stability of Grafted Implant Placement Sites After Sinus Floor Elevation Using a Layering Technique: 10-Year Clinical and Radiographic Results. Int J Oral Maxillofac Implants. 2017;32:1086-1096
- Massimo Del Fabbro, et. al. Systemic Review of Survival Rates for Implants Placed in the Grafted Maxillary Sinus. Int J Periodontics Restorative Dent. 2004;24:565-577
- Nicolaas C. Geurs, et. al. Retrospective Radiographic Analysis of Sinus Graft and Implant Placement Procedures from the Academy of Osseointegration Consensus Conference on Sinus Grafts. Int J Periodontics Restorative Dent. 2001;21:517-523

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
## Reference

Presented by Dr. 손영휘

- Yuki Uchida, et. al. Measurement of Maxillary Sinus Volume Using Computerized Tomographic Images. Int J Oral Maxillofac Implants. 1998;13:811-818
- JH Lee, et. al. Measurement of maxillary sinus volume and available alveolar bone height using computed tomography. Korean J Oral Maxillofac Radiol. 2003;33:35-41
- HW Kim, et al. Measurement of maxillary sinus volume for the placement of graft material – A case control study using CT. J. Kor. Oral Maxillofac. Surg. 2007;33:511-517

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
#### Reference

- Paulo Maló DDS, et. al. All-on-Four™ Immediate-Function Concept with Brånemark System® Implants for Completely Edentulous Mandibles: A Retrospective Clinical Study. Clin Implant Dent Relat Res. 2005 Feb;7 Suppl 1(1, supplement):S88-94
- Sebastian B M Patzelt, et. al. The all-on-four treatment concept: a systematic review. Clin Implant Dent Relat Res. 2014 Dec;16(6):836-55
- David Soto-Penalzo, et al. The all-on-four treatment concept: Systematic review. J Clin Exp Dent. 2017 Mar 1;9(3):e474-e488

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
#### Reference

- GTP(The glossary of prosthodontic terms) ninth edition
- 대한치과보철학회 치과 보철학 용어집 제4판
- Stewarts' clinical removable partial prosthodontics fourth edition
- EB Bae, et. al. A Clinical Retrospective Study of Distal Extension Removable Partial Denture with Implant Surveyed Bridge or Stud Type Attachment. Biomed Res Int. 2017;2017:7140870
- YK Oh, et al. Retrospective clinical evaluation of implant-assisted removable partial dentures combined with implant surveyed prostheses. J Prosthet Dent. 2021 Jul;126(1):76-82
- SH Kang, et. al. Survival rate and clinical evaluation of the implants in implant assisted removable partial dentures: surveyed crown and overdenture. J Adv Prosthodont. 2020 Aug;12(4):239-249

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#### Reference

Presented by Dr. 주현철, 김학후

- Lars Sennerby, 20 years of experience with resonance frequency analysis. Implantologie. 21(1):21-33.
- Vladimir Kokovic, et. al. Immediate vs. early loading of SLA implants in the posterior mandible: 5-year results of randomized controlled clinical trial. Clin Oral Implants Res. 2014 Feb;25(2):e114-9.
- Michael M Bornstein, et. al. Early loading of nonsubmerged titanium implants with a chemically modified sand-blasted and acid-etched surface: 6-month results of a prospective case series study in the posterior mandible focusing on peri-implant crestal bone changes and implant stability quotient (ISQ) values. Clin Implant Dent Relat Res. 2009 Dec;11(4):338-47.
- Serge Baltayan, et. al. The Predictive Value of Resonance Frequency Analysis Measurements in the Surgical Placement and Loading of Endosseous Implants. J Oral Maxillofac Surg. 2016 Jun;74(6):1145-52.
- Pär-Olov Ostman, et. al. Direct implant loading in the edentulous maxilla using a bone density-adapted surgical protocol and primary implant stability criteria for inclusion. Clin Implant Dent Relat Res. 2005;7 Suppl 1:S60-9.
- Daniel Rodrigo, et. al. Diagnosis of implant stability and its impact on implant survival: a prospective case series study. Clin Oral Implants Res. 21, 2010; 255-261.
- L Pagliani, et. al. The relationship between resonance frequency analysis (RFA) and lateral displacement of dental implants: an in vitro study. J Oral Rehabil. 2013 Mar;40(3):221-7.
- Paolo Trisi PhD, et. al. Implant Stability Quotient (ISQ) vs Direct in Vitro Measurement of Primary Stability (Micromotion): Effect of Bone Density and Insertion Torque. J Osteol Biomat. 2010; 1:141-151.
- Stefan Paul Hicklin, et. al. Early loading of titanium dental implants with an intra-operatively conditioned hydrophilic implant surface after 21 days of healing. Clin Oral Implants Res. 2016 Jul;27(7):875-83.
- IH Cho, et. al. A comparative study on the accuracy of the devices for measuring the implant stability. J Adv Prosthodont. 2009 Nov;1(3):124-8.
- R Jaramillo, et. al. Comparative analysis of 2 resonance frequency measurement devices: Osstell Mentor and Osstell ISQ. Implant Dent. 2014 Jun;23(3):351-6.

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#### Reference

Presented by Dr. 박종현

- BS Jeon, et. al. Possibility and Limitation of Applying Occlusal Adjustment on Implant Prosthesis similar to Natural Teeth. Journal of implantology and applied sciences. September 2018. 162-172
- YB Park., Clinical considerations for appropriate occlusion of implant restorations. J Kor DentAssoc. 2015; 53: 418-426
- Shiming Liu, et. al. Influence of occlusal contact and cusp inclination on the biomechanical character of a maxillary premolar: A finite element analysisA review of the literature. J Prosthet Dent 2014;112:1238-1245
- Jae-Hong Lee, et. al. Association between dental implants in the posterior region and traumatic occlusion in the adjacent premolars: a long-term follow-up clinical and radiographic analysis. J Periodontal Implant Sci. 2016 Dec;46(6):396-404
- L A Weinberg. Reduction of implant loading with therapeutic biomechanics. Implant Dent. 998;7(4):277-85
- Okeson JP. Management of temporomandibular disorders and occlusion. 7th ed. St. Louis: Mosby/Elsevier; 2013. p. 443-468
- BS Jeon, et. al. Possibility and Limitation of Applying Occlusal Adjustment on Implant Prosthesis similar to Natural Teeth. Implantology 2018; 22(3): 162-172
- Kaukinen JA, et. al. The influence of occlusal design on simulated masticatory forces transferred to implant-retained prostheses and supporting bone. J Prosthet Dent. 1996 Jul;76(1):50-5.
- The Korean Academy of Stomatognathic Function & Occlusion. Textbook of human jaw function & occlusion. Seoul: Yenang Inc.; 2014 p. 363-7, 387
- R. J. Chapman, et al. Principles of occlusion of implant prostheses: guidelines for position, timing, and force of occlusal contacts. Quintessence Int. 1989; 20: 473-480
- D Wismeijer, et. al. Factors to consider in selecting an occlusal concept for patients with implants in the edentulous mandible. J Prosthet Dent. 1995 Oct;74(4):380-4

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#### Reference

Presented by Dr. 김학후

The logo for OSSTEM IMPLANT features the word "OSSTEM" in a bold, white, sans-serif font with a registered trademark symbol (®) to its upper right. Below it, the word "IMPLANT" is written in a smaller, white, sans-serif font. The background of the entire page is a dark blue-grey color, with a large, abstract graphic on the right side consisting of numerous concentric, overlapping circles in a lighter blue-grey shade, creating a ripple effect.

**OSSTEM<sup>®</sup>**  
**IMPLANT**