

Super convenient 2-step procedure

122 Taper KIT

- Surgical kit optimized for Taper implants
- Reduction of treatment time through increasingly convenient procedures
- Increased user comfort as no cortical drilling is required

Soft Bone

1


Done in One Steps!



Normal Bone

2


Done in Two Steps!



Hard Bone

2

Done in Two Steps!



Excellent Drilling Power



Maximize drilling performance (with outstanding drill bit design and manufacturing technology), skipping intermediate drilling and achieving consistent and fast drilling results.

Preventing Bone Heating

The excellent drilling performance allows intermediate holes to be skipped and this prevents bone heating by maintaining a low drilling temperature.

Preventing Drill Bouncing

The multi-shaped taper drill also acts as a pilot drill, preventing impact during drilling.



Drilling temperature when using

