

# K5 Dental Unit Chair

Our K5 Unit Chair is thoughtfully engineered to elevate both patient experience and clinical workflow. Featuring a slim backrest that minimizes interference with the practitioner's knee space, a comfortable seat designed for longer procedures, and a two-joint headrest operable via knee movement for diverse patient positioning.



## SUPERIOR PATIENT COMFORT

Ergonomic seat and dual-axis headrest accommodate patients of all ages.



## HYGIENIC DELIVERY SYSTEM

Built-in sterilized water supply and pipe cleaning system for sanitary excellence.




## ENHANCED CLINICAL CONTROL


Wireless foot pedal allows hands-free operation, improving workflow efficiency.



## Product Details



**Slim Backrest:** Optimal treatment position by minimizing interference with the user's knee.



**Two-Joint Head Rest:** Designed to treat children, elderly, and patients restricted to wheelchair.



**Sterilizable Doctor's Table-Handle:** Can be separated & sterilized for easy cleaning and contamination.



**Built-in Water Pipe Cleaning System:** Oxygenated water facilitates effective cleaning of the water pipewater pipe.