

N1

The N1 redefines dental imaging convenience. Engineered with a compact, ergonomically balanced “gun” design weighing just 1.61 kg, it allows comfortable one-hand operation - ideal for dynamic dental environments. Despite its portability, it delivers image quality comparable to standard X-ray units (70 kV, 3 mA) with a fine 0.3 mm focal spot. Its high-efficiency battery supports up to 300 images on a single weekly charge, boosting treatment-room flexibility.



EXCEPTIONAL IMAGE QUALITY

Produces diagnostic-grade X-rays equivalent to traditional systems using 70 kV, 3 mA settings and a 0.3 mm focal spot.



ULTRA-LIGHT ONE-HANDED OPERATION

Ergonomic design and ~1.6 kg weight eliminate wrist strain, allowing precise control and comfort.




HIGH-CAPACITY BATTERY

Enables up to 300 shots per charge - recharging only needed once a week with average use.



Product Details



Thanks to 3mA the N1 provides clear, standard quality X-ray quality images within the same radiation time.



Light weight at only 1.61 kg and usable with one hand thanks to the “Gun-Type” design



Long-lasting battery which make the device usable up to 1 week on a single charge (30 images per day)