

TS Implant

The TS (Transcendet Solution) Implant System is one of the most widely used implants among dental professionals worldwide. Designed for optimal primary stability and enhanced osseointegration, this bone-level implant provides predictable outcomes in a wide range of clinical cases.



HIGH PRIMARY STABILITY & EASY INSERTION

Self-tapping, dual-thread design and optimized body shape ensure secure, effortless placement.




FAIR & ACCESSIBLE

High-volume production enables attractive pricing – without compromising on quality.



PROSTHETIC VERSATILITY

Comprehensive prosthetic solutions for both aesthetic and functional restorations.




Product Details




A wide variety of abutments & components maximizes optimal patient care



Aggressive corkscrew threads for faster and stabilized implant insertion



Triple helix cutting edge for a self tapping effect & chip pocket function (cuts & collects bone)



Open Thread helps preventing necrosis on the cervical margin