

# Crest Remover

The Crest Remover is a specialized drill designed to flatten narrow alveolar ridges and establish an accurate entry point for implant osteotomy. Its precise cutting geometry reduces bone irregularities and ensures more predictable guide trajectory.



## ACCURATE RIDGE FLATTENING

Designed to flatten narrow ridges precisely for a more controlled start to implant osteotomy.




## CLEAN AND CONTROLLED CUTTING

Special cutting geometry reduces bone splintering and supports efficient bone removal.




## IMPROVED PROCEDURAL PREDICTABILITY

Creates a smoother crest surface for more stable drill guidance and predictable site preparation.



## Product Details



The blade is designed to minimize bone splintering effectively during the drilling process



Net-shaped blades of the bur speeds up the ridge widening



Creates a smoother crest surface for more stable drill guidance and predictable site preparation