

ENABLES IMPLANT PLACEMENT WITH SIMPLIFIED DRILLING SEQUENCE



122 Taper KIT

Complete implant placement in just **2 simple steps** with the 122 Taper Kit. This user-friendly surgical solution streamlines your workflow, reduces chairside time, and helps **minimize treatment risks** - so you can treat more patients in less time with a **safe, reliable kit**.

» Based on surgical procedure for 4.5 x 10mm implant in soft bone



IMPLANT PLACEMENT IN 2 STEPS

Soft bone implants in 1 step

Normal bone implants in 2 steps

Hard bone implants in 2 Steps



PREVENTS BONE HEATING

Excellent drilling power leads to skipping intermediate drilling which prevents bone heating by maintaining low drilling temperature.



PREVENTING DRILL BOUNCING


Multiple shaped taper drill also functions as pilot drill and this prevents bouncing during drilling.




Surgical Sequence



Soft Bone



Normal Bone



Hard Bone

122 Taper KIT (O122TPK)

Top panel components:

Torque Wrench
TW30B



Depth Gauge
OSDG

