

Assist KIT

Implant surgery often requires multiple instruments for soft tissue access, ridge shaping, bone collection, and implant removal. The Assist Kit **delivers all of these tools in a single solution**, saving time and improving efficiency. With reduced trauma, better precision, and simplified workflows, clinicians can deliver faster, safer, and more predictable implant procedures - while keeping instruments neatly organized.



ALL-IN-ONE CONVENIENCE

Five essential tools in one compact surgical kit.



BONE PRESERVATION & GRAFTING

The Trepine drill's thin 0.4 mm blade minimizes bone loss, while the Autobone collector allows safe autogenous bone harvesting during drilling.



MINIMALLY INVASIVE & PROSTHETIC ACCURACY

The Tissue punch & crest remover reduce healing time and trauma, while the Bone profiler ensures clean seating of abutments and components.



Tools in the Assist Kit



Tissue Punch



Bone Profiler



Trepine Drill



Crestal Remover



Autobone Collector

Assist KIT (OAK)

