

IM-Cure KIT

Designed to combat peri-implant disease, the Osstem IM-Cure Kit equips clinicians with precision instruments to remove biofilm, calculus, and foreign particles from implant surfaces and surrounding tissue. Surgical-grade yet gentle, it helps restore peri-implant health without compromising the implant surface.



SAFE SURFACE CLEANING

Minimizes damage to the implant surface during maintenance of peri-implantitis-affected sites.



DUAL-PROBE DIAGNOSIS

Diagnoses mucositis around the soft tissue without scratching & measures peri-implant pocket depth and extent of peri-implantitis.



VERSATILE SCALERS

Metal SmartScaler: Removes subgingival debris while preserving surface integrity. Plastic Scaler Tip: Safely debrides abutment or crown surfaces without scratching.




Benefits



SmartBrush1: Used for dehiscence & horizontal bone loss



SmartBrush2: Used for vertical bone loss



SmartScaler - Metal: Flexible for easy bending to enhance accessibility



SmartScaler - Plastic: Removal of Peri-Mucositis

IM-Cure KIT (HICK)

