

# MS KIT

The MS Kit is a precision-engineered surgical kit optimized for **implant placement in narrow anterior mandibular ridges** - exceptional for situations with limited interdental bone width. It is specifically designed for seamless MS Implant application in tight anatomical spaces.



## NARROW RIDGE OPTIMIZATION

Tailored for anterior mandibular regions with limited bone width, ensuring accurate and safe MS Implant placement.



## PRECISE DEPTH CONTROL

The included straight drill allows for controlled, accurate drilling depth and reliable implant placement.




## SIMPLIFIED DRILLING WORKFLOW


Featuring an intuitive drilling sequence, the MS Kit streamlines surgical steps, enhancing efficiency and minimizing chair time.



## Examples of Use



Narrow Ridge



Denture



Provisional

# MS KIT (OMSK)

Applicable Products:

MS

Lower panel components:

**Depth Gauge**  
MSDG



**Ratchet Wrench**  
CITQW-1185A

